

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Acid Phosphatase  
**Catalog # :** HA703

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[ACP]
ASSAY CODE	[RATE-A]:[28]-[50]
SAMPLE VOLUME	[20] [5]
R1 VOLUME	[250] [20] [NO]
R2 VOLUME	
WAVELENGTH	[660] [415]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[7200] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Prepare reagents according to package insert instructions.  
Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.  
\*\* Input appropriate bottle size.  
Set K-Factor to 1553

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.  
Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: Albumin

Catalog # : HA901

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[ALB]
ASSAY CODE	[1-POINT]:[8]-[0]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[350] [100] [NO]
R2 VOLUME	
WAVELENGTH	[700] [600]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[1000]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 63

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Alkaline Phosphatase  
**Catalog # :** HA916

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[ALP]
ASSAY CODE	[RATE-A]:[30]-[50]
SAMPLE VOLUME	[6] [3]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[50] [**] [No]
WAVELENGTH	[660] [415]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[11000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 2935. Adjust if necessary

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: ALT

Catalog # : HA926

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[ALT]
ASSAY CODE	[RATE-A]:[30]-[50]
SAMPLE VOLUME	[10] [5]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[50] [50] [No]
WAVELENGTH	[700] [340]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[5500] [DECREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to -5354. Adjust if necessary

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Amylase

**Catalog # :** HA965

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[AMY]
ASSAY CODE	[RATE-A]:[40]-[50]
SAMPLE VOLUME	[10] [2]
R1 VOLUME	[250] [20] [NO]
R2 VOLUME	[50] [20] [No]
WAVELENGTH	[700] [415]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[20000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 3887. Adjust if necessary

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Apolipoprotein A1  
**Catalog # :** A7544

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[APOA1]
ASSAY CODE	[2-POINT]:[23]-[50]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[300] [100] [NO]
R2 VOLUME	[75] [20] [No]
WAVELENGTH	[700] [340]
CALIBRATION	[NONLINEAR] [1] [5]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	[*] [*]
STD (4) CONC-POS	[*] [*]
STD (5) CONC-POS	[*] [*]
STD (6) CONC-POS	[*] [*]
SD LIMIT	[300]
DUPLICATE LIMIT	[500]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagents provided as ready to use Liquids.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.  
Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Apolipoprotein B  
**Catalog # :** A7588

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[APOB]
ASSAY CODE	[2-POINT]:[23]-[50]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[300] [100] [NO]
R2 VOLUME	[75] [20] [No]
WAVELENGTH	[700] [340]
CALIBRATION	[NONLINEAR] [1] [5]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	[*] [*]
STD (4) CONC-POS	[*] [*]
STD (5) CONC-POS	[*] [*]
STD (6) CONC-POS	[*] [*]
SD LIMIT	[500]
DUPLICATE LIMIT	[500]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagents provided as ready to use Liquids.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.  
Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: ASO

Catalog # : A7566

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[ASO]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[15] [15]
R1 VOLUME	[250] [20] [NO]
R2 VOLUME	[50] [20] [No]
WAVELENGTH	[700] [340]
CALIBRATION	[NONLINEAR] [4] [6]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[999]
DUPLICATE LIMIT	[10000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Prepare reagents according to package insert instructions.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.



# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: AST

Catalog # : HA961

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[AST]
ASSAY CODE	[RATE-A]:[30]-[50]
SAMPLE VOLUME	[10] [5]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[50] [50] [No]
WAVELENGTH	[700] [340]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[5500] [DECREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to -5288. Adjust if necessary

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Direct Bilirubin

**Catalog # :** HB936

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[DBIL]
ASSAY CODE	[2-POINT]:[24]-[30]
SAMPLE VOLUME	[6] [3]
R1 VOLUME	[250] [50] [NO]
R2 VOLUME	[65] [20] [No]
WAVELENGTH	[660] [546]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.00] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[32000] [UPPER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Total Bilirubin

**Catalog # :** HB979

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[TBIL]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[4] [2]
R1 VOLUME	[250] [100] [No]
R2 VOLUME	[65] [20] [No]
WAVELENGTH	[600] [546]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[30]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[32000] [UPPER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 1440

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** BUN

**Catalog # :** HB952

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[UN]
ASSAY CODE	[2-POINT]:[28]-[35]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[300] [100] [NO]
R2 VOLUME	[60] [50] [No]
WAVELENGTH	[376] [340]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[7000] [DECREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to -724

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Beta Hydroxybutyrate  
**Catalog # :** H7587

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	BHY
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[8] [4]
R1 VOLUME	[300] [20] [NO]
R2 VOLUME	[50] [20] [No]
WAVELENGTH	[ ] [505]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Calcium (CPC)

**Catalog # :** HC902

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[CA]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[10] [5]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[100] [100] [No]
WAVELENGTH	[700] [600]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[150]
SENSITIVITY LIMIT	[1200]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[32000] [UPPER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 330

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Carbon Dioxide  
**Catalog # :** HC704

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[CO2]
ASSAY CODE	[2-POINT]:[3]-[6]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[300] [20] [NO]
R2 VOLUME	
WAVELENGTH	[376] [340]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [DECREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Mix R1 and R2 together prior to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to -163

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Cholesterol

**Catalog # :** HC910

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[CHOL]
ASSAY CODE	[1-POINT]:[24]-[0]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	
WAVELENGTH	[700] [505]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[1500]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 500

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.



# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** CRP (HS)

**Catalog # :** C7564-40

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[CRP]
ASSAY CODE	[2-POINT]:[28]-[43]
SAMPLE VOLUME	[15] [6]
R1 VOLUME	[150] [20] [NO]
R2 VOLUME	[100] [20] [No]
WAVELENGTH	[800] [570]
CALIBRATION	[NONLINEAR] [4] [6]
STD (1) CONC-POS	[0.00] [1]
STD (2) CONC-POS	[0.05] [*]
STD (3) CONC-POS	[0.15] [*]
STD (4) CONC-POS	[0.50] [*]
STD (5) CONC-POS	[0.75] [*]
STD (6) CONC-POS	[1.0] [*]
SD LIMIT	[999]
DUPLICATE LIMIT	[10000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: CRP (HS) Wide Range

Catalog # : C7568

## CHEMISTRY PARAMETERS:

TEST NAME	[CRP]
ASSAY CODE	[2-POINT]:[28]-[50]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[125] [20] [NO]
R2 VOLUME	[125] [20] [No]
WAVELENGTH	[800] [570]
CALIBRATION	[NONLINEAR] [4] [6]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[5.0] [*]
STD (3) CONC-POS	[20.0] [*]
STD (4) CONC-POS	[40.0] [*]
STD (5) CONC-POS	[160.0] [*]
STD (6) CONC-POS	[320.0] [*]
SD LIMIT	[999]
DUPLICATE LIMIT	[10000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Creatinine

**Catalog # :** HC939

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[CRE]
ASSAY CODE	[2-POINT]:[30]-[50]
SAMPLE VOLUME	[10] [5]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[50] [100] [No]
WAVELENGTH	[570] [505]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[10]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[4500] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 747

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Creatine Kinase  
**Catalog # :** HC922

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[CK]
ASSAY CODE	[RATE-A]:[35]-[50]
SAMPLE VOLUME	[7] [3]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[50] [20] [No]
WAVELENGTH	[376] [340]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[13000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 8306. Adjust if necessary.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Creatine Kinase MB  
**Catalog # :** C7562

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[CK-MB]
ASSAY CODE	[RATE-A]:[35]-[50]
SAMPLE VOLUME	[14] [7]
R1 VOLUME	[300] [100] [NO]
R2 VOLUME	[ ] [20] [NO]
WAVELENGTH	[376] [340]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[13000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 8306. Adjust if necessary.

NON-VALIDATED APPLICATION

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: Fructosamine

Catalog # : F7546

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[FRUC]
ASSAY CODE	[2-POINT]:[45]-[50]
SAMPLE VOLUME	[20] [7]
R1 VOLUME	[250] [*] [NO]
R2 VOLUME	
WAVELENGTH	[700] [546]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.00] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[500]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[32000] [UPPER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Prepare reagents according to package insert instructions.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: G-6-PD

Catalog # : G7583

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	G6PD
ASSAY CODE	[RATE-A]:[15]-[35]
SAMPLE VOLUME	[12] [12]
R1 VOLUME	[350] [20] [No]
R2 VOLUME	
WAVELENGTH	[376] [340]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[999]
DUPLICATE LIMIT	[10000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Prepare reagent by reconstituting 6.0 ml vial with 6.0 mls DH2O. Let dissolve and then add 12 mls R2 reagent to the same vial. This will be the working reagent.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

Prepare samples by mixing 100ul whole blood with 0.9 mls lyse. Let sit 5 minutes. Mix well.

Set K-Factor to 98377

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** GGT

**Catalog # :** HG959

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[GGT]
ASSAY CODE	[RATE-A]:[30]-[35]
SAMPLE VOLUME	[7] [2]
R1 VOLUME	[250] [50] [NO]
R2 VOLUME	[100] [50] [No]
WAVELENGTH	[700] [415]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[5500] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 8140. Adjust if necessary.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.



# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Glucose (Hexokinase)  
**Catalog # :** HG920

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[GLU]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[50] [20] [No]
WAVELENGTH	[376] [340]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[20000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 450

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717  
Test: Hemoglobin A1c  
Catalog # : H7541

## CHEMISTRY PARAMETERS:

TEST NAME	[HbA1c]
ASSAY CODE	[1-POINT]:[50]-[0]
SAMPLE VOLUME	[5] [3]
R1 VOLUME	[180] [20] [NO]
R2 VOLUME	[60] [20] [No]
WAVELENGTH	[ ] [660]
CALIBRATION	[NONLINEAR] [4] [5]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [2]
STD (3) CONC-POS	[*] [3]
STD (4) CONC-POS	[*] [4]
STD (5) CONC-POS	[*] [5]
STD (6) CONC-POS	
SD LIMIT	[999]
DUPLICATE LIMIT	[1000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[-] [-]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

Set K-Factor to 4200

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

Rev. 8-03

Distributed in Mexico:

MYM Laboratory Supply  
California,(619)425-7026

Grupo MOSCARO, S.A. de C.V.  
Mexico, D.F. (Ask MYM)

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Glucose (Oxidase)  
**Catalog # :** HG921

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[GLU]
ASSAY CODE	[1-POINT]:[50]-[0]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[300] [100] [NO]
R2 VOLUME	
WAVELENGTH	[700] [505]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[1100]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 450

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** autoHDL Cholesterol  
**Catalog # :** H7545

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[HDL]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[4] [2]
R1 VOLUME	[300] [**] [NO]
R2 VOLUME	[100] [**] [No]
WAVELENGTH	[700] [600]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[80]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[32000] [UPPER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 600

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Total Iron

**Catalog # :** HI904

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[FE]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[15] [10]
R1 VOLUME	[250] [100] [No]
R2 VOLUME	[50] [100] [No]
WAVELENGTH	[700] [570]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 5160

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: UIBC

Catalog # : 17506

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[UIBC]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[15] [0]
R1 VOLUME	[250] [*] [No]
R2 VOLUME	[50] [*] [No]
WAVELENGTH	[700] [570]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[500] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

R1 = Wkg UIBC Buffer, add 1.7 parts Iron Std (17505-STD) to 10.0 parts UIBC Buffer.

R2 = Iron Color (Cat. No.HI904\_R2)

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 3546

Saline is STD#1(input as 0ug/dl), 500ug/dl Iron standard is STD#2

After calibration record the K-factor. Under the Parameter screen change the "Calibration Type" to K-factor.

Under the Calibration list enter the K-factor as a negative number. You are now ready to run patient samples. The printed result is the UIBC result.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: Beta Hydroxy

Catalog # : H7587

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	BHY]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[8] [4]
R1 VOLUME	[300] [20] [NO]
R2 VOLUME	[50] [20] [No]
WAVELENGTH	[ ] [505]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 256

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.





# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Lactate

**Catalog # :** L7596

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[LAC]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[150] [50] [NO]
R2 VOLUME	[100] [20] [No]
WAVELENGTH	[700] [546]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[1000]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** autoLDL Cholesterol  
**Catalog # :** L7574

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[LDL]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[300] [*] [NO]
R2 VOLUME	[100] [*] [No]
WAVELENGTH	[660] [546]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[1000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 480

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** LDH-L

**Catalog # :** HL956

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[LDH]
ASSAY CODE	[RATE-A]:[35]-[46]
SAMPLE VOLUME	[10] [5]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[100] [20] [No]
WAVELENGTH	[376] [340]
CALIBRATION	[K-FACTOR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[10000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 8864. Adjust if necessary.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

Analyzer: Hitachi 717

Test: Lipase

Catalog # : L7503

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[LIP]
ASSAY CODE	[2-POINT]:[42]-[50]
SAMPLE VOLUME	[5] [3]
R1 VOLUME	[300] [**] [NO]
R2 VOLUME	[100] [**] [No]
WAVELENGTH	[660] [546]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[1000] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Prepare reagents according to package insert instructions.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Magnesium

**Catalog # :** HM929

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[MG]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[4] [2]
R1 VOLUME	[250] [50] [NO]
R2 VOLUME	[250] [50] [No]
WAVELENGTH	[700] [546]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[250]
SENSITIVITY LIMIT	[2000]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 125

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Microalbumin

**Catalog # :** M7562

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[u-ALB]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[10] [10]
R1 VOLUME	[300] [20] [NO]
R2 VOLUME	[100] [20] [No]
WAVELENGTH	[700] [340]
CALIBRATION	[NONLINEAR] [4] [6]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[0.5] [*]
STD (3) CONC-POS	[1.0] [*]
STD (4) CONC-POS	[5.0] [*]
STD (5) CONC-POS	[10.0] [*]
STD (6) CONC-POS	[30.0] [*]
SD LIMIT	[100]
DUPLICATE LIMIT	[10000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 22136

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Microprotein

**Catalog # :** HP782

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[MIPRO]
ASSAY CODE	[1-POINT]:[50]-[0]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[300] [50] [NO]
R2 VOLUME	
WAVELENGTH	[700] [600]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 600

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Phosphorus

**Catalog # :** HP916

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[PHOS]
ASSAY CODE	[2-POINT]:[22]-[37]
SAMPLE VOLUME	[5] [3]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[100] [100] [No]
WAVELENGTH	[700] [340]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[2000]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 170

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.



# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717  
**Test:** Rheumatoid Factor  
**Catalog # :** R7568

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[RF]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[15] [2]
R1 VOLUME	[250] [20] [NO]
R2 VOLUME	[75] [20] [No]
WAVELENGTH	[700] [340]
CALIBRATION	[NONLINEAR] [4] [6]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[999]
DUPLICATE LIMIT	[10000]
SENSITIVITY LIMIT	[0]
ABS LIMIT (INC/DEC)	[32000] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

It is recommended that two levels of control material be assayed daily.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Total Protein

**Catalog # :** HT928

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[TP]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[100] [100] [No]
WAVELENGTH	[700] [570]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[1000]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[-32000] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 475

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Triglycerides

**Catalog # :** HT932

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[TRIG]
ASSAY CODE	[1-POINT]:[50]-[0]
SAMPLE VOLUME	[3] [2]
R1 VOLUME	[300] [100] [NO]
R2 VOLUME	
WAVELENGTH	[700] [505]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[200]
SENSITIVITY LIMIT	[1100]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[0] [LOWER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 875

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.

# Pointe Scientific, Inc.

Instrument Application

**Analyzer:** Hitachi 717

**Test:** Uric Acid

**Catalog # :** HU982

---

---

## CHEMISTRY PARAMETERS:

TEST NAME	[UA]
ASSAY CODE	[2-POINT]:[24]-[50]
SAMPLE VOLUME	[7] [3]
R1 VOLUME	[250] [100] [NO]
R2 VOLUME	[50] [20] [No]
WAVELENGTH	[660] [546]
CALIBRATION	[LINEAR] [0] [0]
STD (1) CONC-POS	[0.0] [1]
STD (2) CONC-POS	[*] [*]
STD (3) CONC-POS	
STD (4) CONC-POS	
STD (5) CONC-POS	
STD (6) CONC-POS	
SD LIMIT	[0.1]
DUPLICATE LIMIT	[100]
SENSITIVITY LIMIT	[700]
ABS LIMIT (INC/DEC)	[0] [INCREASE]
PROZONE LIMIT	[32000] [UPPER]
EXPECTED VALUE	[*] [*]
PANIC VALUE	[*] [*]
INSTRUMENT FACTOR	[1.0]

### REAGENT PREPARATION

Reagent is ready to use.

Please consult the package insert for complete instructions on the use of this reagent.

### PROCEDURAL NOTES

\* Indicates user defined parameter.

\*\* Input appropriate bottle size.

Set K-Factor to 568

It is recommended that two levels of control material be assayed daily. Reorder PSI Chemistry Controls Cat.# C7590-50 & C7591-50.

Hitachi 717™ is a Registered Trademark of Nissei Sangyo Co. Ltd., Japan.