

Installation Instructions for Abbott Cell-Dyn® 1600 & 1700 Reagent Set

Required ONLY when changing to CDS reagents from other non-CDS products.

For proper performance DO NOT mix CDS reagents with non-CDS products

To achieve optimal performance with CDS NextGeneration™ Reagents, it is important to completely prime the Abbott Cell-Dyn® instrument. To do so, simply follow the steps below:

To prime Diluent and Detergent

1. Remove pickup tubes from diluent and detergent containers, press REAGENT PRIME
2. When reagent empty alarm is activated, place pickup tubes in appropriate CDS Reagents (Diluent—CDS 501-067 or 501-069, Detergent—501-064 or 501-066), then press CLEAR ALARM to continue reagent prime
3. Press REAGENT PRIME three (3) more times

To prime Lyse

1. Remove lyse pickup line from container, press REAGENT PRIME
2. When reagent empty alarm is activated, place pickup tube in CDS Lytic Agent (CDS 501-080 or CDS 501-077), then press CLEAR ALARM to continue reagent prime
3. Press Lyse prime three (3) more times
4. Run three (3) BACKGROUND CYCLES to complete reagent exchange

The instrument is now completely primed with

CDS NextGeneration™ Reagents.

Proceed with normal instrument operation, complete reagent log.

If you experience any problems or have any questions, please call 800-453-3328.