

# Rotores para Centrifuga C5 USA

---

Rotor de 4 cubetas con aditamento azul para insertar 6 tubos por cubeta.

Gira 24 tubos de 3 – 10ml. con muestras de sangre.

Gira 8 tubes de 15 ml. para orina

Max 4,000 rpm / 2800g



## 4-Place Bucket Rotor

BLUE 6-place test tube inserts

-Spins 24 3ml-10ml blood test tubes

-Spins Crit-Carriers

-Spins 8 15ml conical urine tubes

**Rotor de 8 plazas para tubos de 15ml.  
Max 5,000 rpm / 4250g**



**8-Place 15ml Rotor**  
Spins 8 3ml-15ml test tubes  
Spins Crit Carriers  
Max 5,000rpm / 4250g

**Rotor de 4 cubetas con dispositivo negro para insertar tubos de 50ml. para orina y semen.**



- BLACK 4-place 50ml inserts**
- Spins 4 50ml tubes for urine and semen
  - Spins 4 Arthrex ACP double syringe
  - Spins StemLogix stem cell tubes
  - Spins StemLogix PRP tubes
  - Spins VetStem 30ml PRP tubes

**Rotor de 4 plazas para tubos de 50ml.  
de orina, semen y PRP  
Max. 5000rpm / 4250g.**



**4-Place 50ml Rotor**

Spins 50ml tubes for urine, semen, and PRP

Spins StemLogix stem cell tubes

Spins StemLogix PRP tubes

Spins VetStem 30ml PRP tubes

Max 5,000rpm / 4250g